

The Symmetry Medical MultiPak® is a modular tray system designed to organize and protect your delicate instruments. Standard MultiPak tray systems are available in Radel® R, aluminum and stainless steel.

Our highly versatile design offers a variety of nesting inserts and accessories to create customized configurations utilizing off the shelf components. Standard accessories include a variety of silicone mats, silicone bar patterns, dividers, posts and wire racks.

Our engineers can further customize your MultiPak with customized colors and materials to further tailor the MultiPak to accommodate virtually any instrument set.

MultiPak is also compatible with the majority of sterilization methods currently available.





5000-01-475P shown



5000-50-000P shown




5000-50-175P shown



5023-03-E shown



5023-50-E shown



5010-03-180P shown

MultiPak® Base Units and Lids

- Large Base Unit Part #: 5000-01-475P
9.2 x 20 x 4.6 inches
(234 x 508 x 117mm)

- Large Base Unit Part #: 5000-01-375P
9.2 x 20 x 3.6 inches
(234 x 508 x 91mm)

- Large Base Unit Part #: 5000-01-275P
9.2 x 20 x 2.6 inches
(234 x 508 x 66mm)

- Lid for 5000 Base Part #: 5000-50-000P
9.7 x 20.5 x 0.2 inches
(246 x 520 x 5mm)

- Elevated Lid for Base Part #: 5000-50-175P
9.7 x 20 x 1.4 inches
(246 x 508 x 36mm)

MultiPak® Insert Trays and Lids

- 2/3 Insert Tray for 5000 Base Part #: 5023-03-E
8.7 x 12.8 x 1.3 inches
(221 x 325 x 33mm)

- 2/3 Clear Lid for 5023-03-E Insert Tray Part #: 5023-50-E
8.8 x 13.1 x 0.2 inches
(224 x 333 x 5mm)

- Full Insert Tray for 5000 Base Part #: 5010-03-180P
8.7 x 19.5 x 1.7 inches
(221 x 495 x 43mm)

- Lid for 5010 Insert Tray Part #: 5010-50-000P
9.1 x 19.9 x 0.2 inches
(231 x 506 x 5mm)

- Elevated Lid for 5010 Insert Tray Part #: 5010-50-175P
9.1 x 19.9 x 0.2 inches
(231 x 506 x 5mm)

- Long Insert Tray for 5000 Base Part #: 5012-03-180P
4 x 19.5 x 1.7 inches
(101.6 x 495 x 43mm)

- Lid for 5012 Insert Tray Part #: 5012-50-000P
4.3 x 19.9 x 0.2 inches
(109 x 506 x 5mm)



5122-03-375P shown



5022-50-175P shown



5023-50-000P shown



5006-01-390P shown



5006-50-000P shown

MultiPak® Insert Trays and Lids

- Insert Tray for 5000 Base Part #: 5022-03-180P

8.7 x 9.5 x 1.7 inches
(221 x 241 x 43mm)

- Lid for 5022 Insert Tray Part #: 5022-50-000P

9 x 9.8 x 0.2 inches
(229 x 249 x 5mm)

- Deep Insert Tray for 5000 Part #: 5122-03-375P

8.7 x 9.5 x 3.3 inches
(221 x 241 x 84mm)

- Elevated Lid for 5022/5122 Insert Tray Part #: 5022-50-175P

8.4 x 9.1 x 1.4 inches
(213 x 231 x 36mm)

- Insert Tray for 5000 Base Part #: 5023-03-180P

8.7 x 12.7 x 1.7 inches
(221 x 323 x 43mm)

- Lid for 5023 Insert Tray Part #: 5023-50-000P

9 x 13 x 0.2 inches
(229 x 330 x 5mm)

- Deep Insert Tray for 5000 Base Part #: 5123-03-375P

8.7 x 12.7 x 3.3 inches
(221 x 353 x 84mm)

- Elevated Lid for 5023/5123 Insert Tray Part #: 5023-50-175P

8.4 x 12.4 x 1.4 inches
(213 x 315 x 36mm)

MultiPak® Jr.

- Base Part #: 5006-01-390P

5.5 x 12.8 x 3.9 inches
(140 x 325 x 98mm)

- Insert for 5006 Base Part #: 5006-03-150P

5.1 x 11.2 x 1 inches
(130 x 285 x 25mm)

- Lid for 5006 Base Part #: 5006-50-000P

5.8 x 12.9 x 0.2 inches
(147 x 328 x 5mm)

The MultiPak® Endoscopic Kit is a highly versatile tray system designed to organize and protect your valuable endoscopic instruments.

The modular design, coupled with a variety of MultiPak components, permits the custom configuration of virtually any instrument set.

MultiPak® Endoscopic Kit

• MultiPak 20" Endoscopic Kit Part #: 5500


Kit includes: 5000-01-475P Base, 5000-50-175P Lid, 5111 Wire Rack, 5022 Insert Tray, (4) 5200 Bar Holders, (4) 5233 "T" Brackets, and Hardware.*

* Silicone parts not included. Reference MultiPak Accessories for complete list of component options.



MultiPak® Silicone Small Diamond Grid Mats

- Mat for 5000 Base Part #: 5000-04SD
8.9 x 19.7 inches
(226 x 500mm)
- Mat for 5000 Base w/5110 Rack Part #: 5110-04SD
8.9 x 18 inches
(226 x 457mm)
- Mat for 7026 Base Part #: 5023-04SD
9.4 x 25 inches
(239 x 635mm)
- Mat for 5010 Insert Tray Part #: 5010-04SD
8.4 x 19.2 inches
(213 x 488mm)
- Mat for 5012 Insert Tray Part #: 5012-04SD
3.7 x 19.2 inches
(94 x 489mm)
- Mat for 5022 Insert Tray Part #: 5022-04SD
8.4 x 9.1 inches
(213 x 231mm)
- Mat for 5023 Insert Tray Part #: 5023-04SD
8.4 x 12.5 inches
(213 x 318mm)



MultiPak® Silicone Large Diamond Grid Mats

- Mat for 5000 Base Part #: 5000-04D
8.9 x 19.7 inches
(226 x 500mm)
- Mat for 5000 Base w/5110 Rack Part #: 5110-04D
8.9 x 19 inches
(226 x 457mm)
- Mat for 5006 Base Part #: 5006-04D
5.4 x 12.5 inches
(137 x 318mm)
- Mat for 5010 Insert Tray Part #: 5010-04D
8.4 x 19.2 inches
(213 x 488mm)
- Mat for 5012 Insert Tray Part #: 5012-04D
3.7 x 19.2 inches
(94 x 489mm)
- Mat for 5022 Insert Tray Part #: 5022-04D
8.4 x 9.1 inches
(213 x 231mm)
- Mat for 5022 Insert Tray & 5108 Rack Part #: 5108-04D
8.3 x 9 inches
(211 x 229mm)
- Mat for 5023 Insert Tray Part #: 5023-04D
8.4 x 12.5 inches
(213 x 318mm)
- Mat for 5023 Insert Tray & 5109 Rack Part #: 5109-04D
8.3 x 12.4 inches
(211 x 315mm)
- Mat for 5023-E Insert Tray & 5109 Rack Part #: 5023-04E-D
8.4 x 12.5 inches
(213 x 318mm)

MultiPak® Silicone Fine Pin Mats

- Mat for 5000 Base Part #: 5000-04
8.9 x 19.7 inches
(226 x 500mm)

- Mat for 5000 Base w/5110 Rack Part #: 5110-04
8.9 x 18 inches
(226 x 457mm)

- Mat for 5006 Base Part #: 5006-04
5.4 x 12.5 inches
(137 x 318mm)

- Mat for 7024 Base Part #: 7024-04
9.4 x 23 inches
(239 x 584mm)

- Mat for 5010 Insert Tray Part #: 5010-04
8.4 x 19.2 inches
(213 x 488mm)

- Mat for 5012 Insert Tray Part #: 5012-04
3.7 x 19.2 inches
(94 x 489mm)

- Mat for 5022 Insert Tray Part #: 5022-04
8.4 x 9.1 inches
(213 x 231mm)

- Mat for 5023 Insert Tray Part #: 5023-04
8.4 x 12.5 inches
(213 x 318mm)

- Mat for 5023-03-E Insert Tray Part #: 5023-04-E
8.4 x 12.5 inches
(213 x 318mm)



5110 shown with 5200, sold separately

MultiPak® Instrument Racks

- Wire Rack Instrument Holder for 5000 Base Part #: 5110
8.7 x 19.7 x 3.3 inches
(222 x 500 x 83mm)

- Wire Rack Instrument Holder for 5000 Base Part #: 5111
(Requires Elevated Lid Part #: 5000-50-175P)
8.7 x 19.7 x 3.8 inches
(222 x 500 x 95mm)

- Wire Rack with Handle, Instrument Holder for 5000 Base Part #: 5112
(Requires Elevated Lid Part #: 5000-50-175P)
8.7 x 19.7 x 3.9 inches
(222 x 500 x 100mm)


- Wire Rack with Handle, Instrument Holder for 5000 Base Part #: 5113
(Requires Elevated Lid Part #: 5000-50-175P)
8.7 x 19.7 x 4.3 inches
(222 x 500 x 109mm)

- Wire Rack with Handle for 7024 Base Part #: 7113
8.7 x 24.4 x 4.4 inches
(222 x 618 x 112mm)

- Wire Rack with Handle for 7026 Base Part #: 7513
9.34 x 24.5 x 5.5 inches
(237 x 622 x 140mm)

- Wire Rack for 5022 & 5122 Insert Tray Part #: 5108
7.9 x 8.9 x 3 inches
(200 x 227 x 77mm)

- Wire Rack for 5023 & 5123 Insert Tray Part #: 5109
7.9 x 12.9 x 3.2 inches
(200 x 312 x 82mm)



5113 shown

MultiPak® Cases















MultiPak® Accessories*

- Celcon Nuts & Bolts (qty: 10) Part #: 5120 
- Short Divider (qty: 2) Part #: 5235
Compatible with 5000 Base Only, Requires 5215
9.1 x 2.8 inches 
- Long Divider (qty: 2) Part #: 5236
Compatible with 5000 Base Only, Requires 5215
19.9 x 2.8 inches 
- Holder, Slotted for 5235 & 5236 Only Part #: 5215
1.5 inches 
- Short Divider Part #: 5237
For use with 5010 Insert Trays 
- Long Divider Part #: 5238
For use with 5010 Insert Trays 
- Endoscopic Silicone Bar Part #: 5200N
For use with Wire Racks 5108 and 5109
6.5 inches 
- Endoscopic Silicone Bar Part #: 7200
For use with Wire Rack 7513
7 inches 
- Endoscopic Silicone Bar Part #: 5200
For use with Wire Racks 5100 - 5113, 7113
7 inches 
- Holder, Slotted Plastic for Silicone Bars #: 5210
2.9 inches 
- Holder, Slotted Plastic for Silicone Bars #: 5211
7 inches 
- Metal Carrying Handle Kit #: 5130
For use with 5000-01, 5010-03 and 5022-03 
- Metal Screw, Black (qty: 50) Part #: 5121
For use with 5231, 5233
0.25 inches 
- Metal Screw, Black (qty: 50) Part #: 5122
For use with 5241 Post
0.375 inches 

- Nameplate (blank) #: 5216
4.7 x 2.5 inches 
- Plastic Holding Post Part #: 5241
1.41 inches 
- Short Post, Movable Part #: 5250
1.24 inches 
- Long Post, Movable Part #: 5251
2 inches 
- "L" Shelf Bracket Part #: 5231) 
- Stainless Steel "L" Shelf Bracket Part #: 455231
Compatible only with MultiPak III 
- "T" Shelf Bracket Part #: 5233 
- Blank Bar Part #: 5311
2.9 inches 
- Blank Bar (priced per foot) Part #: 5310
12 inches 
- 3 Slot Instrument Bar 3-8mmØ #: 5330
2.9 inches 
- 10 Slot Instrument Bar 3-8mmØ #: 5370
7 inches 
- 3 Slot Instrument Bar 5mmØ #: 5320
2.9 inches 
- 3 Slot Instrument Bar 8-15mmØ #: 5331
2.9 inches 
- 8 Slot Instrument Bar 8-15mmØ #: 5371
7 inches 
- 3 Slot Instrument Bar 10mmØ #: 5334
2.9 inches 
- 8 Slot Instrument Bar 10mmØ #: 5374
7 inches 

*MultiPak Silicone Accessories include plastic holding bar and screws; excluding 5310


MultiPak® Silicone Accessories*

- Double 3 Slot Instrument Bar, 5 & 10mmØ Part #: 5332
2.9 inches 
- Double 9 Slot Instrument Bar, 5 & 10mmØ Part #: 5372
7 inches 
- 4 Slot Instrument Bar Part #: 5312
2.9 inches 
- 5 Slot Instrument Bar, 3mmØ Part #: 5335
2.9 inches 
- 12 Slot Instrument Bar, 3mmØ Part #: 5375
7 inches 
- 5 Slot Instrument Bar, 5mmØ Part #: 5333
2.9 inches 
- 12 Slot Instrument Bar, 5mmØ Part #: 5373
7 inches 
- 6 Slot Instrument Bar, 3mmØ Part #: 5300
7 inches 
- 6 Slot Instrument Bar, 6mmØ Part #: 5301
7 inches 
- 2 Slot Scope Holder, 6mmØ Part #: 5343
Open 15mmØ, Closed 12mmØ
2.9 inches 
- Double Ended Scope Holder Part #: 5341
Open 12mmØ
2.9 inches 
- Double Ended Scope Holder Part #: 5340
Open 12mmØ, Closed 8mmØ
2.9 inches 
- Scope Holder, Open End 12mmØ Part #: 5353
1.5 inches 
- Scope Holder, Closed End 12mmØ Part #: 5352
1.5 inches 

*MultiPak Silicone Accessories include plastic holding bar and screws; excluding 5310

MultiPak® Silicone Small Diamond Grid Mats

- Mat for 5000 Base Part #: 5000-04SD
8.9 x 19.7 inches
(226 x 500mm)
- Mat for 5000 Base w/5110 Rack Part #: 5110-04SD
8.9 x 18 inches
(226 x 457mm)
- Mat for 7026 Base Part #: 5023-04SD
9.4 x 25 inches
(239 x 635mm)
- Mat for 5010 Insert Tray Part #: 5010-04SD
8.4 x 19.2 inches
(213 x 488mm)
- Mat for 5012 Insert Tray Part #: 5012-04SD
3.7 x 19.2 inches
(94 x 489mm)
- Mat for 5022 Insert Tray Part #: 5022-04SD
8.4 x 9.1 inches
(213 x 231mm)
- Mat for 5023 Insert Tray Part #: 5023-04SD
8.4 x 12.5 inches
(213 x 318mm)



MultiPak® Silicone Large Diamond Grid Mats

- Mat for 5000 Base Part #: 5000-04D
8.9 x 19.7 inches
(226 x 500mm)
- Mat for 5000 Base w/5110 Rack Part #: 5110-04D
8.9 x 19 inches
(226 x 457mm)
- Mat for 5006 Base Part #: 5006-04D
5.4 x 12.5 inches
(137 x 318mm)
- Mat for 5010 Insert Tray Part #: 5010-04D
8.4 x 19.2 inches
(213 x 488mm)
- Mat for 5012 Insert Tray Part #: 5012-04D
3.7 x 19.2 inches
(94 x 489mm)
- Mat for 5022 Insert Tray Part #: 5022-04D
8.4 x 9.1 inches
(213 x 231mm)
- Mat for 5022 Insert Tray & 5108 Rack Part #: 5108-04D
8.3 x 9 inches
(211 x 229mm)
- Mat for 5023 Insert Tray Part #: 5023-04D
8.4 x 12.5 inches
(213 x 318mm)
- Mat for 5023 Insert Tray & 5109 Rack Part #: 5109-04D
8.3 x 12.4 inches
(211 x 315mm)
- Mat for 5023-E Insert Tray & 5109 Rack Part #: 5023-04E-D
8.4 x 12.5 inches
(213 x 318mm)