









Urology from HSW

Cystoscopes 2,7 mm

Description 	HD	Part number	Working length	Total length
Cystoscope 2,7 mm autoclavable 0°	✓	7199202700	270 mm	340 mm
Cystoscope 2,7 mm autoclavable 30°	✓	7199202730	270 mm	340 mm

Cystoscopes 4,0 mm

Description 	HD	Part number	Working length	Total length
Cystoscope 4,0 mm autoclavable 0°		8300003265	300 mm	355 mm
Cystoscope 4,0 mm autoclavable 12°		8300003267	300 mm	355 mm
Cystoscope 4,0 mm autoclavable 30°		8300003270	300 mm	355 mm
Cystoscope 4,0 mm autoclavable 70°		8300003271	300 mm	355 mm