

All-in-one Med-cam system

Ultra-slim Video Nasolaryngoscope

High-quality HD Camera Head

All-purpose system for ENT applications



All-in-one – All you need

This all-in-one Med-cam system has been developed in close cooperation with clinical users in order to provide maximum benefits for both patients and users. The ultra-slim Video Nasolaryngoscope provides an extremely gentle insertion, meaning minimal discomfort and time in procedure. The high-quality HD camera head offers all possibilities and enables optimal image quality.

CMOS Video Nasolaryngoscope – Optimal cost / performance ratio



Integrated LED
for optimal operating expense



Ergonomic design
for effortless use



Ultra-slim-line
for gentle insertion




Precise motion
for perfect control

CMOS HD Camera Head – All possibilities

Ergonomic design
for effortless use

All-purpose system
for ENT applications

HD sensor
for optimal image quality



| Video Nasolaryngoscope with integrated LED illumination | |
|---|----------------|
| Article no. | 830 001 1825 |
| Working length | approx. 300 mm |
| Diameter distal end | 3.2 mm |
| Diameter insertion tube | 3.4 mm |
| Direction of view | 0° |
| Field of view | 90° |
| Bending range of distal tip | 130° up / down |

| Video Nasolaryngoscope, light guide version | |
|---|--------------|
| Article no. | 830 001 2731 |

| Camera Head | |
|----------------------|-------------------------------|
| Article no. | 830 001 3118 |
| Coupler | |
| Article no. | 830 001 8854 |
| Med-cam control unit | |
| Article no. | 830 001 2748 |
| Input voltage range | 100 - 240 V AC, 50 / 60 Hz |
| Output | 1x USB for image storage, DVI |