

## *All-in-one Med-cam system*

*Ultra-slim Video Nasolaryngoscope*

*High-quality HD Camera Head*

*All-purpose system for ENT applications*



# All-in-one – All you need

This all-in-one Med-cam system has been developed in close cooperation with clinical users in order to provide maximum benefits for both patients and users. The ultra-slim Video Nasolaryngoscope provides an extremely gentle insertion, meaning minimal discomfort and time in procedure. The high-quality HD camera head offers all possibilities and enables optimal image quality.

## CMOS Video Nasolaryngoscope – Optimal cost / performance ratio



**Integrated LED**  
for optimal operating expense



**Ergonomic design**  
for effortless use



**Ultra-slim-line**  
for gentle insertion




**Precise motion**  
for perfect control

## CMOS HD Camera Head – All possibilities

**Ergonomic design**  
for effortless use

**All-purpose system**  
for ENT applications

**HD sensor**  
for optimal image quality



Video Nasolaryngoscope with integrated LED illumination	
Article no.	830 001 1825
Working length	approx. 300 mm
Diameter distal end	3.2 mm
Diameter insertion tube	3.4 mm
Direction of view	0°
Field of view	90°
Bending range of distal tip	130° up /down

Video Nasolaryngoscope, light guide version	
Article no.	830 001 2731

Camera Head	
Article no.	830 001 3118
Coupler	
Article no.	830 001 8854
Med-cam control unit	
Article no.	830 001 2748
Input voltage range	100 - 240 V AC, 50/60 Hz
Output	1x USB for image storage, DVI